



BEFOUR

PROFESSIONAL HEALTHCARE SCALES

Since 1979, Befour, Inc. has been designing & manufacturing precision digital scales for the most demanding environments. Befour's line of Professional Healthcare Scales have been specifically engineered with both the care provider and patient in mind. Every scale has an unsurpassed combination of advanced features, safety and ease of use.

EXPECT MORE.

BEFOUR

PROFESSIONAL HEALTHCARE SCALES

...the important information you need - in a comprehensive guide

Befour's line of professional healthcare scales will serve just about any need within a healthcare setting. Using the same display console across the product line offers ease of operation and requires minimal training for both medical staff and bio-med departments.




PORTABLE PLATFORM, COLUMN & HANDRAIL SCALES



Model	MX115	MX160	MX170	MX305	MX306	MX308	MX377	MX310
Weight Capacity	550 lbs	750 lbs	1000 lbs	750 lbs	750 lbs	750 lbs	750 lbs	1000 lbs
Platform Size (w x d)	13" x 13"	16" x 18"	24" x 20"	16" x 18"	16" x 18"	16" x 18"	16" x 18"	22" x 25"
System Weight	12 lbs	14 lbs	25 lbs	17 lbs	31 lbs	31 lbs	44 lbs	46 lbs
Resolution	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb
Accuracy	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb
Portability	Carry	Carry	Carry Handle		Tilt & Roll	Tilt & Roll		Tilt & Roll
Handrail					Included	Included	Included	Included
Integrated Height Rod								
Alkaline Battery Power	●	●	●	●	●	●	●	●
Optional AC Adapter	●	●	●	●	●	●	●	●
Body Mass Index (BMI)	●	●	●	●	●	●	●	●
Stand-On	●	●	●	●	●	●	●	●

MEASUREMENT STATIONS



Model	MX805 / MX905	MX807	MX808 / MX908	MX810 / MX910	MX877 / MX977
Weight Capacity	750 lbs	750 lbs	750 lbs	1000 lbs	750 lbs
Platform Size (w x d)	16" x 18"	16" x 18"	16" x 18"	22" x 25"	16" x 18"
System Weight	17 lbs	25 lbs	33 lbs	46 lbs	44 lbs
Resolution	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb
Accuracy	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb
Portability		Tilt & Roll	Tilt & Roll	Tilt & Roll	
Handrail		Included	Included	Included	Included
Integrated Height Rod	●	●	●	●	●
 Opti-Height Optical Stadiometer System	●		●	●	●
Alkaline Battery Power	●	●	●	●	●
Optional AC Adapter	●	●	●	●	●
Body Mass Index (BMI)	●	●	●	●	●
Stand-On	●	●	●	●	●

THE SCALE PRODUCT LINE COMPARISON GUIDE




All scales feature



Technology for accurate weight readings

Industry Leading - 3 year warranty

WHEELCHAIR SCALES



Model	MX420	MX450	MX450D	MX480	MX480D	MX490D
Weight Capacity	750 lbs	1000 lbs	1000 lbs	1000 lbs	1000 lbs	1000 lbs
Platform Size (w x d)	Variable	28" x 32"	28" x 32"	32" x 34"	32" x 34"	32" x 41"
Ramps	Integral	Single	Dual	Single	Dual	Dual
System Weight	24 lbs	65 lbs	65 lbs	87 lbs	87 lbs	100 lbs
Resolution	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb	0.1 lb
Accuracy	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb	± 0.1 lb
Portability	Carry	Tilt & Roll	Tilt & Roll	Fold & Roll	Fold & Roll	Fold & Roll
Handrail		Included	Included	Included	Included	Included
Alkaline Battery Power	●	●	●	●	●	●
Optional AC Adapter	●	●	●	●	●	●
Body Mass Index (BMI)	●	●	●	●	●	●
Stand-On		●	●	●	●	●
Chair (optional use)	●	●	●	●	●	●
Wheelchair	●	●	●	●	●	●

INFANT & SPECIALTY SCALES



Model	MX308CHR	MX700	MX202	MX282
Weight Capacity	500 lbs	1000 lbs	45 lbs	60 lbs
Platform Size (w x d)	21" x 17" seat	32" x 50"	24" x 16-1/2"	39" x 18"
System Weight	64 lbs	120 lbs	27 lbs	37 lbs
Resolution	0.1 lb	0.1 lb	1 gram	0.02 lb
Accuracy	± 0.1 lb	± 0.1 lb	± 1 gram	± 0.02 lb
Portability	Tilt & Roll		Handles	
Handrail	Included	Included		
Alkaline Battery Power	●	●	●	●
Optional AC Adapter	●	●	●	●
Body Mass Index (BMI)	●	●		●
Stand-On	Converts to Stand on	●		
Chair	●	●		
Wheelchair		●		

Why choose a Befour Precision Digital Scale?

- ✓ Proudly Made in the USA
- ✓ Industry-leading 3 year warranty
- ✓ Quality construction for repetitive use in demanding environments
- ✓ Active handrails & integrated digital height rods included in many models
- ✓ MotionLock Technology engineered into every scale
- ✓ Flexible, customer-centric manufacturer with a passion for service and quality
- ✓ And EVERYTHING listed below!



Many scales, but only ONE display console

Did you know that Befour's line of Professional Healthcare Scales all use the SAME 7-Button Display? Streamline your scales and seamlessly move from model to model. Without the need to retrain your Biomed teams and caregivers, you'll save both TIME & MONEY!

Every Display Console Features ...

Smart-View:	Show the information you want, weight only or weight, height & BMI
Programmable User Settings:	Tailor the display to your needs including auto-off time, "Hands Free" activation & more!
LB/KG Lockout:	Toggle back and forth or LOCK-IN a desired unit of measure
Body Mass Index:	Display quickly and automatically calculates BMI
ECO Weigh Battery System:	4 D-cell alkaline batteries (included) provide 100,000+ weighings
Battery-Life Indicator:	3-segment indicator provide plenty of notice before battery replacement is needed
Push-Button Calibration:	Scales arrive 100% calibrated but if ever needed, recalibration is just a few simple steps
EMR/EHR Ready:	Connectivity via USB serial output or optional wireless plug-in module
Stored Tare:	Store the weight of a frequently used item such as a wheelchair ... And more!

BEFOUR

Befour Inc. | 102 Progress Drive, Saukville, WI 53080 | 800.367.7109 | f: 262.284.5966 | befour.com | sales@befour.com