

BEFOUR

PROFESSIONAL HEALTHCARE SCALES

THE Exam Room Scale

MODEL MX377

Befour's **PREMIER** physician scale - THE Exam Room Scale - offers every feature you'd want in a stationary patient scale.

Stability, safety and technology are at the forefront of engineering. An innovative 'dual axis tip-proof' design and integrated "active" handrail - which comes standard - provide unmatched safety. This free-standing, rock-solid, scale incorporates all-metal construction to deliver robust quality in a sleek, space-saving package.

The base model (MX377) of THE Exam Room Scale features a 750 pound capacity, Befour's proprietary MotionLock technology, and ECO-Weigh power management system providing more than 100,000 weight readings on a single set of batteries. All Befour scales are **built to last**, backed by an 'industry leading' **3-year warranty** and are **Made in the USA**.

EMR READY



HANDS-FREE
WEIGH

FEATURES & BENEFITS

Integrated "Active" Handrail: A STANDARD feature, allowing patients to grasp handrail for support while weight is accurately calculated

'Dual Axis Tip-Proof' Design: Unique, free-standing design eliminates all front-to-back and side-to-side tipping for unmatched stability and patient security

Display Console: LB/KG lockout, programmable settings, BMI output, auto zero, recall, reweigh ability, and more



MotionLock automatically facilitates quick & accurate weight, even for the most unsteady user – no intervention of the care provider or patient required



DISPLAY CONSOLE

Display is used across all Befour's medical scales: Seamlessly move from one model to another with ease

Displays weight, height, and BMI on a single screen

EMR connectivity via serial output or other options

Programmable User Settings can be set to auto-off, "Hands-Free" Weighing, etc.

'ECO Weigh' battery system provides 100,000 weighings per set of (4) D-cell alkaline batteries (included)

LB/KG Lockout: Toggle between or set to display only your desired unit of measure

3-segment battery life indicator

Large 1" LCD digits

No AC outlet or recharging needed

THE Exam Room Scale

MODELS MX377



SCALE SPECIFICATIONS

SCALE FOOTPRINT	16" W x 22" D x 46" H
PLATFORM DIMENSIONS	16" W x 18" D x 2" H
SYSTEM WEIGHT	41-44 lbs
WEIGHT CAPACITY	750 lbs (340 kg)
RESOLUTION	0.1 lb (0.05 kg)
ACCURACY	± 0.1 lb (0.05 kg)
POWER SOURCE	4 D-cell alkaline batteries (incl) *optional AC adapter available (P/N: 03049-08)
BATTERY LIFE	100,000 weighings
CONNECTIVITY	USB serial or other options
DISPLAY READOUT	1" LCD digits
WARRANTY	3 years

SHIPPING SPECIFICATIONS

PACKAGED WEIGHT	~ 45-48 lbs
PACKAGE DIMENSIONS	~ Box (48 x 20" x 10")
SHIPPING METHOD	UPS Package

DIMENSIONS



DUAL AXIS TIP-PROOF DESIGN *

Befour's exclusive 'dual axis tip-proof' engineering offers rock-solid stability, in a free-standing design, producing an incredibly secure scale.

* Design reduces chance of tipping under normal usage.



BEFOUR

Befour Inc. | 102 Progress Drive, Saukville, WI 53080 | 800.367.7109 | f: 262.284.5966 | befour.com | sales@befour.com