

BEFOUR

SPORTS & FITNESS SCALES

Portable Battery-Powered Scale with BMI

Model PS-6700

Large platform and battery-operated display make this a favored choice! This portable battery-operated scale delivers a quick and accurate weight. The PS-6700 scale is extremely durable yet surprisingly lightweight.

Designed with functional basics in mind, this ruggedly built scale with MotionLock technology provides quick, reliable and extremely accurate weight measurements, even for unsteady users.

As the top-selling national brand in team sports, Befour weigh-in scales combine accuracy, reliability and durability in a truly portable design.

A clear choice for athletic trainers, home health care professionals, weight-loss clinics, and sporting events such as wrestling, martial arts, football and more. No outlet or re-charging needed!

FEATURES & BENEFITS:

Pin point accuracy of ± 0.1 lb.

High Capacity: 500 lbs.

Fully assembled: Arrives ready to use; simply insert the batteries (included)

Display Console Set-up: Quick insert/detach connections

Push Button Calibration: Set calibration via keypad entry without the need of tools

Rugged and durable all-metal construction

MOTIONLOCK Accurately calculates the weight of an unsteady user

Plug & Weigh Exclusive cable system can be easily replaced if damaged



DISPLAY CONSOLE



- **Displays weight, height and BMI** on a single screen
- **Programmable User Settings:** can be set to auto-off (12 seconds to 25 minutes), always-on, etc.
- **LB/KG Lockout:** Change setup to display only your desired unit of measure
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed (batteries last 100,000+ weighings)
- Optional AC Adapter available: P/N 03049-06
- Auto zero, recall, and reweigh

Optional Carry Cases:



SC-1816 Soft Carry Case

Extra pockets for documents & supplies; includes shoulder strap



HC-1825 Hard-Shell Carry Case

(lockable) Heavy-duty polyethylene hard-shell case with custom-cut foam interior for added safety during transport.

SPECIFICATIONS:

Platform Dimensions.....	16" x 18" x 1-3/4"
System Weight	15 lbs
Weight Capacity.....	500 lbs (225 kg)
Resolution.....	0.1lb (0.05 kg)
Accuracy.....	± 0.1 lb (0.05kg)
Power Source.....	4 D-cell alkaline batteries
Battery Life.....	100,000+ weighings
Display Readout.....	1" LCD digits
Warranty.....	3 years



Befour, Inc.

102 Progress Drive · Saukville, WI 53080 · P: 800-367-7109 · F: 262-284-5966
Befour.com · sales@befour.com