

BEFOUR

SPORTS & FITNESS SCALES

Model WH-1080

Calculate weight, height and BMI with Befour's transportable WH-1080 scale. WH-1080 is a truly versatile scale that is ideal for most weighing needs. Standard features include sturdy, integrated, wide handrail and Opti-Height Optical Stadiometer System. This exclusive feature enables fast and accurate height measurements through state-of-the-art design and expert engineering. MotionLock technology providing stability for unsteady users and pin-point accuracy of height and weight measurement.

Transport is easy...simply tilt the scale back on its heavy-duty wheels and roll to the next location.

Battery powered—no outlet or recharging needed!

FEATURES & BENEFITS:

Pin point accuracy of $\pm 0.1\text{lb}$

High capacity: 650 lbs

Body Mass Index: Automatically calculates and displays with BMI

Tilt & Roll: Designed for ultimate portability

Push Button Calibration: Set calibration via keypad entry without the need of tools

Rugged and durable all-metal construction

Display Console: LB/KG lockout, programmable settings, BMI output, and more highlighted in right-hand column

Optical Height Measurement: Simply lower the foldable headpiece to rest on top of head, press HEIGHT & measurement occurs instantaneously. Height can be viewed in selectable units and displays automatically



Accurately calculate the weight of an unsteady user



DISPLAY CONSOLE FEATURES



- **Displays weight, height and BMI** on a single screen
- **Programmable User Settings:** can be set to auto-off (12 seconds to 25 minutes), always-on, etc.
- **LB/KG Lockout:** Change setup to display only your desired unit of measure
- 3-segment battery life indicator
- Large 1" LCD digits
- No AC outlet or recharging needed (batteries last 100,000 weighings)
- Optional AC Adapter available: PN 03049-06
- Auto zero, recall, and reweigh



Powered by:

Opti-Height

Optical Stadiometer System

Height rod is both digital and integrated into scale for automatic BMI calculation. Unique design has extended measuring range with no exposed cables and rotates 360° for easy positioning of headpiece

SPECIFICATIONS:

Height Rod Range.....	25-94" (63-238cm)
Platform Dimensions.....	16" x 18" x 1 3/4"
System Weight.....	33 lbs
Resolutions.....	0.1 lb (0.05 kg)
Accuracy.....	± 0.1 lb (0.05 kg)
Power Source.....	4 D-cell Alkaline Batteries
Battery Life.....	100,000 weighings
Display Readout.....	1" LCD digits
Warranty.....	3 years



Befour, Inc.

102 Progress Drive · Saukville, WI 53080 · P: 800-367-7109 · F: 262-284-5966
Befour.com · sales@befour.com